

BT6II *back-to-basics*

Orbital Shaker with
Rubber Mat Tray.



7677-10115
7740-01010

Reciprocating Shaker
with Rubber Mat Tray.



7677-30115
7677-31010

1 Position Stirrer



7767-40115



7677-50115

4 Position Stirrer



7677-10115
7740-02525

Orbital Shaker with
Stainless Steel,
Sticky Strip Tray.



7677-20115
7677-21000

Rocker with Optional
Stacking Tray.

BTBII – BACK-TO-BASICS Rocker Platform

The BTBII Rocker Platform provides an economical, easy to use solution for light duty applications. It is designed for general laboratory products and procedures including gel staining, blot washing and filter assays using flasks, beakers, slides, tubes, etc. It includes a 10x12" platform tray (unit is not compatible with Bellco's standard interchangeable trays). An optional stacking tray is available to expand the capacity of the rocker without sacrificing bench space.

Specifications:

Speed Range: 8-40 rpm

Tilt Angle: +/- 12.5°

Load Capacity: 10 lbs. max – evenly distributed.

Electrical Rating: 115V+/-10 VAC @60Hz line voltage
230V+/-10 VAC @50/60Hz line voltage



Stock

7677-20115

7677-20230

7677-21000

Description

Rocker Platform (115V)

Rocker Platform (230V)

Optional Stacking Tray
(3 1/2" space between trays)

BTBII – BACK-TO-BASIC ORBITAL SHAKER

The BTBII Orbital Shaker provides an economical, easy to use solution for light duty applications. It is designed for general laboratory products and procedures requiring orbital movement; including hybridization and gel agitation. Trays are sold separately.

Specifications:

Speed Range: 30-250 rpm

Orbit: 1" (2.5cm)

Load Capacity: 10 lbs. max – evenly distributed

Electrical Rating: 115V+/-10 VAC @ 60 Hz line voltage
230V+/-10 VAC @50/60 Hz line voltage



Stock

7677-10115

7677-10230

Description

BTBII-Orbital Shaker – 115V

BTBII-Orbital Shaker – 230V

Recommended Trays:

Stock

7740-01010

7740-01212

7740-02525

7740-40000

Description

10"x10" Platform tray (w/mat)

12"x12" Platform tray (w/mat)

10"X10" SST Tray (to be used w/sticky tape)

Sticky Strip Tape

BTBII – BACK-TO-BASICS RECIPROCATING SHAKER

The Reciprocating Shaker provides an economical, easy to use solution for light duty applications. It is designed for general laboratory products and procedures requiring reciprocating stroke movement; including hybridization and gel agitation. Trays are sold separately.

Specifications

Speed Range: 10-200 oscillations per minute

Stroke: 1"

Load Capacity: 10 lbs. max – evenly distributed

Electrical Rating: 115V+/-10 VAC @ 60 Hz line voltages
230V+/-10 VAC @ 50/60Hz line voltage



Stock #	Description
7677-30115*	BTBII-Reciprocating Shaker – 115V
7677-30230*	BTBII-Reciprocating Shaker – 230V

* Consists of: Base unit with 24V in-line power supply (See below for Trays)

Recommended Trays for Reciprocating Units Only:

Stock #	Description
7677-31010	10"x10" Reciprocating Platform tray (w/mat)
7677-32010	10"x10" SST Reciprocating Tray
7740-40000*	Sticky Strip Tape

* Sticky Strip Tape is for SST Platform Only (7677-32010)

Vortex-Genie® 1

Specifications

Supplied with: Pop-Off Cup

Dimensions (base): 165 x 122 x 165mm (D x W x H)
(6.5 x 4.8 x 6.5in)

Weight: 4 Kg (8.8 lb)

Warranty: 2-year manufacturer warranty



Details

- More Violent Vortexing Action than any other vortex mixer.
- Simple to use. Touch and go. Violent high-speed ONLY mixing action eliminates operator variability.
- Affordable. Vortex-Genie quality at an even more economical price.
- Aggressive. Attacks problem applications like no other vortexer can. Vortexing ability has been enhanced
- for maximum high-speed touch mixing efficiency. Reduce vortexing time!
- Solid metal casting and stainless steel construction. Does not "walk".
- Beakers, Flasks or Multiple Tube Vortexing with optional 3" Platform Attachment.

Stock #	Description
	Vortex-Genie 1 Touch Mixer, 120V, 60Hz 0.65 amps
	Vortex-Genie 1 Touch Mixer, 230V, 50Hz 0.5 amps – No Plug
	Vortex-Genie 1 Touch Mixer, 230V, 50Hz 0.5 amps – European Plug
	Vortex-Genie 1 Touch Mixer, 230V, 50Hz 0.5 amps – British Plug
	Vortex-Genie 1 Touch Mixer, 230V, 50Hz 0.5 amps – Swiss Plug

BTBII – BACK-TO-BASICS MAGNETIC STIRRER, 1 & 4 Positions

The BTBII Magnetic Stirrer provides an economical, easy to use solution for light duty applications. It is designed for both cell culture in spinner flasks as well as general laboratory products and procedures requiring the use of magnetic stirring equipment with minimal heat transfer.

Specifications

Speed Range: 50-550 rpm

Load: 1 Position: up to a 500ml flask
4 Position: either (2) 500ml flasks or (4) 250 ml flasks

Electrical Rating: 115V+/-10 VAC @60Hz line voltage
230V+/-10 VAC @50/60Hz line voltage



Stock #	Description
7677-40115*	1 Position 115V
7677-40230*	1 Position 230V
7677-50115*	4 Position 115V
7677-50230*	4 Position 230V

*Consists of: Base unit with 24V in-line power supply



BELLCO GLASS, INC.



Contact: (856) 691-1075 · (800) 257-7043 · FAX: (856) 691-3247
Bellco Glass, Inc. · 340 Edrudo Road · Vineland, NJ · 08360
E-Mail: cservice@bellcoglass.com · Internet: <http://www.bellcoglass.com>